



neo ice®

EQUINE FORMULA RECHARGE SOLUTION

1200mL

Menthol-free Ohne Menthol Sans menthol



COOL & COMPRESS

Competition Safe and effective

Cools without refrigeration

Draws heat for 2½ hours & can reduce limb temperatures by 10°C when used as directed.

Use to recharge depleted Neo-ice Equine Formula Ice-wraps. Apply to horses leg after strenuous activity to reduce inflammation & treat acute injuries. Works with immediate effect.

New formula Compliant with the 2010 guidelines of
Fédération Equestre Internationale (FEI) & Deutsche
Reiterliche Vereinigung, Fédération Equestre Nationale (FN).



INSTRUCTIONS

How to use to recharge Neo-ice Equine Formula Ice-wraps:

1. Place dry Ice-wrap into a plastic bag or vessel.
2. Remove the outer and inner caps above the bottle measuring cylinder and squeeze 80mL into cylinder.
3. Pour Solution onto depleted Ice-wrap by tipping bottle in a clock-wise fashion, and replace bottle caps.
4. Squeeze Ice-wrap for 30-60 seconds to evenly distribute Recharge Solution.

How to apply and use Neo-ice Equine Formula Ice-wrap:


1. Dampen with Neo-ice Equine Formula Recharge Solution just prior to application.
2. Wrap & fasten Ice-wrap around the horses leg.
3. Leave exposed to air - do not cover.
4. For best results, keep the horse immobile and in a well-ventilated environment.

Neo-ice Ice-wraps can be left in place until dry. Wash & dry between applications. Keep bottle tightly sealed and in a cool place when not in use.

CAUTION: Do not use continuously for more than 24 hours: Keep out of the reach of children. Do not apply to broken or damaged skin as contained alcohol will cause stinging. If skin irritation occurs, wash skin and apply moisturizing cream. If splashed in eyes irrigate thoroughly with water and if irritation remains seek medical advice. Do not drink. If swallowed seek medical advice and show label to the Doctor. If swellings do not subside seek veterinary advice.

Contains: Water, Ethanol (denat.), Menthyl lactate, PEG-40 hydrogenated castor oil, CI42090 (blue dye), Triethanolamine.

www.neo-ice.com



Euco Technology Co. Ltd.
www.eucotech.com
Made in Taiwan



5060106031013

Thermal image demonstrates
Neo-ice effectiveness.